
# OctoMag

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
## System Components

### MagnebotiX OctoMag magnetic field and gradient generation unit

- Maximum working volume diameter: 5 cm
- Max. magnetic field strength (center): 50 mT in any direction
- Max. magnetic field gradient (center): 2 T/m
- Dimensions (field generator) W x D x H: 560 mm x 560 mm x 390 mm
- Weight (field generator unit) 125 kg
- Dimensions (outer frame) W x D x H: 640 mm x 940 mm x 1350 mm
- total weight: 180 kg
- inclination angle of field generation unit: 0°, 30° (shown), 45°
- The OctoMag magnetic field generator has no moving parts and can be used over a wide range of temperatures

### MagnebotiX extended power unit (ECB-820-ex)

- Number of channels: 8
- Max. current per channel: 18 A
- Max. voltage per channel: 96 V
- Max. magnetic field frequency: 10 Hz
- Max. power consumption: 2016 W
- External power requirements: 110 - 240 VAC / 50 - 60 Hz
- Dimensions (W x D x H) : 450 mm x 330 mm x 175 mm
- Weight: 12 kg
- Connection to control computer: Ethernet

The control computer and MBX pro software are included with the system.

System safety features: over-temperature and over-current monitoring of all coils; controller shutdown upon loss of communication with coil unit; external emergency shutdown switch.

The following are compatible with the system but not included

- Microscope or optics are chosen according to the experiments to be performed. Dissecting microscopes or long working distance zoom tubes are frequently used. A basic mounting platform is provided.
- Basler USB3 cameras are compatible with our MBX pro software (plug-and-play). Tested products include Basler USB3 ACE or Pulse series. Users can interface their own cameras through the open ROS environment if desired
- A computer monitor is not supplied with the system. We recommend use of a 24" monitor.
- SpaceNavigator (3Dconnexion), PlayStation3 controller (input is provided to the MBX pro control software for open loop control).

The principle of the magnetic field generating technology is protected by patents WO 2011029592 A1 and WO 2013127516 A1 to ETH Zurich and sub-licensed to MagnebotiX AG. For more details on the OctoMag systems see OctoMag: An Electromagnetic System for 5-DOF Wireless Micromanipulation, M. P. Kummer, J. J. Abbott, B. E. Kratochvil, R. Borer, A. Sengul and B. J. Nelson. IEEE Transactions on Robotics, **26** (6), 1006 – 1017, 2010.

The MBX pro control software is based on software originally developed in the Institute of Robotics and Intelligent Systems, Multiscale Robotics Lab, ETH Zürich ("daedalus") and has been adapted to the Robot Operating System (ROS).

*All specifications are subject to change without notice.*