




### Content of the slide

Each Argo-HM slide contains 16 fluorescent patterns.





# Target PAT-AG03-EM1-A1 This pattern consists of concentric circles with increasing radii from 25 µm to 325 µm with a step of 25 µm, plus an extra circle with a radius of 12.5 µm, featuring a target.



1µm step
PAT-AG03-EM1-I1
This pattern consists of twice 11
empty cylinders embedded at
different depths, like two crossing
stairs, surrounded by four pillars.

3D Crossing stairs

The step is: 1 µm.




Field of rings

PAT-AG03-EM1-B1

This pattern consists of a matrix of 45×45 rings, separated by 15 µm, on a total field of 660 µm × 660 µm. The field of rings is surrounded by eight landmarks and exhibits an 8 µm long cross in its center.



3D Crossing stairs 0.75 µm step
PAT-AG03-EM1-I2
This pattern consists of twice 11 empty cylinders embedded at different depths, like two crossing stairs, surrounded by four pillars. The step is: 0.75 µm.



4×4 Intensity gradation
PAT-AG03-EM1-C1
This pattern consists of sixteen
8.5 μm-wide squares having
different fluorescence intensity
levels following a linear evolution,
organized in a 4×4 matrix.



3D Crossing stairs 0.5 µm step PAT-AG03-EM1-I3

This pattern consists of twice 11 empty cylinders embedded at different depths, like two crossing stairs, surrounded by four pillars. The step is: 0.5 µm.



2×16 Intensity Gradation

PAT-AG03-EM1-D1

This pattern consists of twice sixteen 22.5 µm × 1.5 µm rectangles having different fluorescence intensity levels following a linear evolution, organized in a 2×16 matrix.



3D Crossing stairs 0.25 µm step

PAT-AG03-EM1-I4

This pattern consists of twice 11 empty cylinders embedded at different depths, like two crossing stairs, surrounded by four pillars. The step is: 0.25 µm.

上海昊量光电设备有限公司 地址:上海市徐汇区远中产业园三期6号楼3层

电话:400-6888-532 网站: www.auniontech.com





Gradually spaced lines descending

PAT-AG03-EM1-E3
This pattern consists of pairs of 50 µm-long lines whose spacing gradually increases, from 100 nm to 700 nm, with a step of 50 nm.
One set of lines is present:



Repositionning crosses
PAT-AG03-EM1-H1
The repositioning crosses are 40
µm long.




Gradually spaced lines horizontal

PAT-AG03-EM1-E1
This pattern consists of pairs of 50 µm-long lines whose spacing gradually increases, from 100 nm to 700 nm, with a step of 50 nm.
One set of lines is present: horizontal.




Sphere
PAT-AG03-EM1-G1
This pattern consists of three circles with a diameter of 50 µm in different orthogonal planes, featuring the equator and two meridians of a sphere.



# Gradually spaced lines ascending

PAT-AG03-EM1-E4
This pattern consists of pairs of 50 µm-long lines whose spacing gradually increases, from 100 nm to 700 nm, with a step of 50 nm.
One set of lines is present: ascending (+ 45°).




### Matrix of crosses

PAT-AG03-EM1-F1 This pattern consists of a matrix of  $6\times6$  crosses, having a length of 5  $\mu$ m and a step of 8  $\mu$ m, surrounded by a 50  $\mu$ m-wide frame.

The crosses are composed of vertical lines that are in the same plane, and by horizontal lines, going gradually deeper within the glass.

The spacing between the vertical and horizontal lines gradually increases, from 0.1  $\mu$ m to 3.6  $\mu$ m, with a step of 0.1  $\mu$ m.



# Gradually spaced lines vertical

PAT-AG03-EM1-E2 This pattern consists of pairs of 50 µm-long lines whose spacing

50 µm-long lines whose spacing gradually increases, from 100 nm to 700 nm, with a step of 50 nm. One set of lines is present: vertical.



### Word ARGOLIGHT

电话:400-6888-532

网站: www. auniontech. com

PAT-AG03-EM1-J1
This pattern consists of the letters forming the company name "Argolight", and surrounded by a 220 µm × 50 µm frame.

Lifetime warrantied fluorescence presence

Dimensions:75x25x1.5 mm

Materials: Anodized aluminum enclosure with an AG03 glass core

Excitation range:continuum 250-650 nm

Emission range:continuum from the excitation wavelength plus 15 nm,to 800 nm

Immersion medium compatibility:dry, oil: no limitation,water objectives: less than 20 min at a time

Storage conditions:room temperature (10-40  $^{\circ}$ C)and under normal relative humidity(20-70  $^{\circ}$ RH) Imaging technology compatibility:any fluorescence-based imaging except depletion-based technology and multiphoton imaging

Light exposure damage threshold:50 GW/cm2 irradiance (peak or average)



## **Analyze**

More than 12 automated quality tests, several tens of relevant metrics.



FIELD UNIFORMITY



FIELD DISTORTION



LATERAL CO-REGISTRATION ACCURACY



LINE SPREAD FUNCTION



RING SPREAD FUNCTION



LATERAL RESOLUTION



OPTICAL SECTIONING STRENGTH



STAGE REPOSITIONING REPEATABILITY



STAGE DRIFT DURING TIMELAPSE



STAGE DRIFT DURING Z STACKING



ACCURACY OF 3D RECONSTRUCTION

电话:400-6888-532

网站: www. auniontech. com



INTENSITY RESPONSE



SPECTRAL RESPONSE



POWER METER



POINT SPREAD FUNCTION



ACCURACY OF CO-REGISTRATION



UNIFORMITY OF FIELD