

## ISTEQ Broadband plasma light sources

- XWS-65: Basic configuration
- XWS-Dual port: Two light outputs
- XWS-R: High power light source
- XWS-X: Ultra high brightness
- XWS-30: Compact light source
- Hyperchromator: Tunable light source

# XWS-65 broadband plasma light source

## XWS-65 laser pumped plasma ultrabright broadband light source

ISTEQ's XWS-65 light source product has been specially developed to be used for a variety of applications, including spectroscopy, high resolution microscopy, thin - film measurement, surface metrology and others. This source is based on cutting edge technology, covered by EU and US patents.




### Application fields:

- Absorption and fluorescence spectroscopy
- Diagnostics systems in microelectronics - contamination and defect control
- Surface metrology, ellipsometry and scatterometry
- Microscopy, including confocal and fluorescence
- Optical component testing
- Detectors in chromatography, microfluidics, lab-on-a-chip, droplet spectrometers, cytofluorimeters, etc

### Main advantages:

- CW laser plasma discharge
- Broad spectral range: 190 – 2500 nm
- High spectral brightness: up to 50 mW/(mm<sup>2</sup>·sr·nm)
- High temporal and spatial stability: STD < 0.15%
- Long life time due to electrodeless operation: 10,000 hours
- The small dimensions of the emitting volume considerably expand the range of XWS applications
- External source control and parameters monitoring via Software, Windows GUI

## Spectral brightness of XWS-65 light source in UV and VIS spectral region



# XWS-65 broadband plasma light source

## XWS-65 performance

Spectral range	190 - 2500nm for UV configuration, 250 - 2500nm for OFR configuration
Spectral brightness	Up to 50mW/(mm <sup>2</sup> ·sr·nm)
Output power	Up to 3 W free space Up to 0.5 W via fiber
Lamp medium	Xenon
Emitting body size	250×500um
Lifetime	10,000 hours
Temporal and spatial stability	STD < 0.15%

## Optical design

Output NA by default	0.4, up to 0.55 upon request
External optic interface by default	C-mount
Optional output interface	Thorlabs SM1, 30mm cage and more
Fiber interface (only for FCU version)	SMA or FC

## Optional configurations

Source spectrum	UV or Ozone free
Light output	Free space or fiber coupled
Optical head cooling unit	Air or water cooling
Power Supply Unit (PSU) cooling system	Air or water cooling

## Additional

Connection to PC/Laptop	Ethernet (Web interface), COM-port (RS232)
Interlock	Db-15 connector
Remote plasma control	Db-15 connector

## System dimensions and weight


Optical head Free Space	130 × 110 × 74mm, 1.3kg
Power supply unit	351 × 172 × 232mm, 8kg

## Facility requirements


Electrical	100-240V, 50/60Hz
Gas purging (only for UV configuration)	Nitrogen or Argon purging, 1l/min

# XWS-65 broadband plasma light source

## Optical head (Free Space) dimensions:




## Power Supply Unit (PSU) dimensions:






# XWS-65 FCU broadband plasma light source


## XWS-65 OFR FCU configuration:



## XWS-65 UV FCU configuration:



# XWS-Dual port: Free space and FCU versions



**XWS Dual port Free Space**

**XWS Dual port FCU**


**Special features:**

- Dual beam light source with identical parameters for each port
- Developed on XWS-65 base
- Spectral brightness: up to 50 mW/(mm<sup>2</sup>·sr·nm)
- High stability: STD <0.15%

**FCU version:**

- SMA or FC fiber interface
- Output power up to 0.5W after each fiber

## Spectral brightness of XWS-Dual port light source in UV and VIS spectral region



# XWS-Dual port: Free space and FCU versions

## XWS-Dual port performance

Spectral range	190 - 2500nm for UV configuration, 250 - 2500nm for OFR configuration
Spectral brightness	Up to 50mW/(mm <sup>2</sup> ·sr·nm)
Output power	Up to 3 W free space per port Up to 0.5 W via fiber per port
Lamp medium	Xenon
Emitting body size	250×500um
Lifetime	10,000 hours
Temporal and spatial stability	STD < 0.15%

## Optical design

Output NA by default	0.4, up to 0.55 upon request
External optic interface by default	C-mount
Optional output interface	Thorlabs SM1, 30mm cage and more
Fiber interface (only for FCU version)	SMA or FC

## Optional configurations

Source spectrum	UV or Ozone free
Light output	Free space or fiber coupled
Optical head cooling unit	Air cooling
Power Supply Unit (PSU) cooling system	Air or water cooling

## Additional

Connection to PC/Laptop	Ethernet (Web interface), COM-port (RS232)
Interlock	Db-15 connector
Remote plasma control	Db-15 connector

## System dimensions and weight


Optical head FCU	130 × 154 × 74mm, 2kg
Power supply unit	351 × 172 × 232mm, 8kg

## Facility requirements


Electrical	100-240V, 50/60Hz
Gas purging (only for UV configuration)	Nitrogen or Argon purging, 1l/min

# XWS-Dual port FCU

## Optical head dimensions:




## Power Supply Unit (PSU) dimensions:






# XWS-Dual port Free Space


## Optical head dimensions:



## Power Supply Unit (PSU) dimensions:




## XWS-R high power light source



### Special features:

- High power version: almost doubled output power/brightness in comparison to default XWS-65
- Spectral brightness: up to 90 mW/(mm<sup>2</sup>·sr·nm)
- Can be done in both: Free space and FCU configuration
- Up to 1W of output power after a fiber for FCU version

### Spectral brightness of XWS-R light source in UV and VIS spectral region



# XWS-R high power light source

## XWS-R performance

Spectral range	190 - 2500nm for UV configuration, 250 - 2500nm for OFR configuration
Spectral brightness	Up to 90mW/(mm <sup>2</sup> ·sr·nm)
Output power	Up to 5 W free space per port Up to 1 W via fiber per port
Lamp medium	Xenon
Emitting body size	250×700um
Lifetime	10,000 hours
Temporal and spatial stability	STD < 0.15%

## Optical design

Output NA by default	0.4, up to 0.55 upon request
External optic interface by default	C-mount
Optional output interface	Thorlabs SM1, 30mm cage and more
Fiber interface (only for FCU version)	SMA or FC

## Optional configurations

Source spectrum	UV or Ozone free
Light output	Free space or fiber coupled
Optical head cooling unit	Air cooling
Power Supply Unit (PSU) cooling system	Air or water cooling

## Additional

Connection to PC/Laptop	Ethernet (Web interface), COM-port (RS232)
Interlock	Db-15 connector
Remote plasma control	Db-15 connector

## System dimensions and weight


Optical head FCU	130 × 130 × 74mm, 2kg
Power supply unit	351 × 172 × 232mm, 8kg

## Facility requirements


Electrical	100-240V, 50/60Hz
Gas purging (only for UV configuration)	Nitrogen or Argon purging, 1l/min

# XWS-R high power light source

## Optical head dimensions:



## Power Supply Unit (PSU) dimensions:





## XWS-X high UV light source


ISTEQ's XWS-X light source product has been specially developed for those customers, who needs a powerful light source with a high UV generation.

### Special features:

- Significantly higher brightness across the whole spectral range in comparison to XWS-65
- Maximum spectral brightness up to 120mW/(mm<sup>2</sup>·sr·nm)
- X6 brightness improvement at UV range below 250nm in comparison to XWS-65
- Output configuration: Free space or FCU
- Available in Dual Port configuration
- Available in XR (High power) configuration with a maximum brightness up to 200mW/(mm<sup>2</sup>·sr·nm)



### Spectral brightness of XWS-X light source in UV and VIS spectral region



# XWS-X high power light source

## XWS-X performance

Spectral range	190 - 2500nm for UV configuration, 250 - 2500nm for OFR configuration
Spectral brightness	Up to 120mW/(mm <sup>2</sup> ·sr·nm)
Output power	Up to 3 W free space Up to 0.6 W via fiber
Lamp medium	Xenon
Emitting body size	250×500um
Lifetime	10,000 hours
Temporal and spatial stability	STD < 0.15%

## Optical design

Output NA by default	0.4, up to 0.55 upon request
External optic interface by default	C-mount
Optional output interface	Thorlabs SM1, 30mm cage and more
Fiber interface (only for FCU version)	SMA or FC

## Optional configurations

Source spectrum	UV or Ozone free
Light output	Free space or fiber coupled
Optical head cooling unit	Air cooling
Power Supply Unit (PSU) cooling system	Air or water cooling

## Additional

Connection to PC/Laptop	Ethernet (Web interface), COM-port (RS232)
Interlock	Db-15 connector
Remote plasma control	Db-15 connector

## System dimensions and weight


Optical head FCU	130 × 110 × 74mm, 1.8kg
Power supply unit	351 × 172 × 232mm, 8kg

## Facility requirements


Electrical	100-240V, 50/60Hz
Gas purging (only for UV configuration)	Nitrogen or Argon purging, 1l/min

# XWS-X broadband plasma light source

## Optical head (Free Space) dimensions:



## Power Supply Unit (PSU) dimensions:



## XWS-30 compact light source




ISTEQ's XWS-30 light source product has been specially developed for those customers, who needs super compact broadband light source with low heat dissipation, keeping the plasma brightness high.

### Special features:

- Unique concept of a compact **"all in one"** source
- Very compact: 110x110x120mm, no external chiller
- Spectral brightness: up to 40 mW/(mm<sup>2</sup>·sr·nm)
- Output configuration: Free space or FCU
- Full system control by Laptop/PC via USB-RS485 adapter

### Spectral brightness of XWS-30 light source in UV and VIS spectral region





# XWS-30 high power light source

## XWS-30 performance

Spectral range	190 - 2500nm for UV configuration, 250 - 2500nm for OFR configuration
Spectral brightness	Up to 40mW/(mm <sup>2</sup> ·sr·nm)
Output power	Up to 1.5 W free space Up to 0.4 W via fiber
Lamp medium	Xenon
Emitting body size	100×200um
Lifetime	10,000 hours
Temporal and spatial stability	STD < 0.15%

## Optical design

Output NA by default	0.5, up to 0.55 upon request
External optic interface by default	C-mount
Optional output interface:	Thorlabs SM1, 30mm cage and more
Fiber interface (only for FCU version)	SMA or FC

## Optional configurations

Source spectrum	UV or Ozone free
Light output	Free space or fiber coupled
Optical head cooling unit	Air cooling

## Additional

Connection to PC/Laptop	COM-port (RS485) or USB
Interlock	Lemo FG


## System dimensions and weight

Optical head	110 × 110 × 120mm, 2kg
--------------	------------------------

## Facility requirements

Electrical	100-240V, 50/60Hz
Gas purging (only for UV configuration)	Nitrogen or Argon purging, 1l/min

# XWS-30 compact light source



## XWS specifications

	XWS-30	XWS-65	XWS-R	XWS-X
Spectral range for UV/OFR:	190 - 2500nm / 250 - 2500nm			
Maximum spectral brightness	40mW/(mm <sup>2</sup> ·nm·sr)	50mW/(mm <sup>2</sup> ·nm·sr)	90mW/(mm <sup>2</sup> ·nm·sr)	120mW/(mm <sup>2</sup> ·nm·sr)
Output power	Up to 1.5W free space Up to 0.4W via fiber	Up to 3W free space Up to 0.5W via fiber	Up to 5W free space Up to 1W via fiber	Up to 3W free space Up to 0.6W via fiber
Emitting body size	100×200um	250×500um	250×700um	250×500um
Lifetime	10,000 hours			
Temporal and spatial stability	STD < 0.15%			
<b>Optical design</b>				
Output NA (by default)	0.5	0.4	0.4	0.4
Maximum NA (upon request)	0.55			
External optic interface	C-mount, 30 Thorlabs cage			
<b>Optional configurations</b>				
Source spectrum	UV or Ozone free			
Light output	Free space or fiber coupled			
<b>Additional features</b>				
External control	COM-port (RS-485)	COM-port (RS-232)		
Interlock / Distant plasma control	Lemo FGG	Db-15 connector		
<b>Dimensions</b>				
Optical head dimensions	138x110x148mm	130x75x106mm	130x110x106mm	130x75x106mm
Power supply dimensions	NA (No power supply)	351x175x232mm	351x175x232mm	351x175x232mm

## Available configurations

	Output configuration		Cooling system			
	Dual Port Output	UV/OFR FCU	Water Cooled Optical Head	Air Cooled PSU	Water Cooled PSU	Air Cooled PSU
XWS-30	✘	✔	✘	✔	N/A	N/A
XWS-65	✔	✔	✔	✔	✔	✔
XWS-R	✘	✔	✔	✔	✔	✔
XWS-X	✔	✔	✔	✔	✔	✔
XWS-XR	✘	✔	✔	✔	✔	✔

## Custom design

Upon a request that is possible to modify ISTEQ serial product to meet customer requirements.

上海昊量光电设备有限公司 Phone: 4006-888-532 WeChat: Auniontech Website: www.auniontech.com E-mail: info@auniontech.com



# Hyperchromator – tunable light source




The Hyperchromator - «Mountain Instruments» monochromator designed by for ISTEQ's laser-pumped plasma light source XWS-30. The Hyperchromator optimally uses the advantages of the pointed plasma light source by using a mirror with an aperture up to f/1.5 to collect the light and ensure maximum light throughput without a limiting entrance slit. It is optimized for applications that require a fast and tunable point light source.

## Special features:

- Different configurations available
- Fast optics, up to f/1.5 for highest throughput
- Homogenous output distribution due to a proprietary design
- Etendue-matched to ISTEQ XWS-30
- Broad tunable range from DUV to NIR
- No input slit
- Built-In Shutter
- Easy to use Software, Windows GUI
- LabView integration upon request

## Power Hyperchromator II 400 $\mu$ m fiber with NA=0.22





# Hyperchromator – tunable light source

## Hyperchromator performance

Optical input	ISTEQ XWS-30 light source, directly coupled (optionally other sources)
Optical output	Fused silica fiber, SMA or FC, 100-600 $\mu\text{m}$ core or free beam output with adjustable slit or various collimator options. Spectral power monitoring on request.
Wavelength range	185 – 2500 nm*
Aperture	f/1.5 or f/2 (depending on required resolution and light output)
Bandwidth	1-10 nm FWHM*
Output power	Up to 800 $\mu\text{W}$ (grating at blaze wavelength, 6 nm bandwidth and 400 $\mu\text{m}$ fiber)
Reproducibility	0.1 nm
Scanning speed	40-100 nm/s*
Control interface	USB/RS-232, LabVIEW™-based GUI, various external control options
Dimensions and Weight	47 x 45 x 25cm (WxDxH); 16kg

\*: depends on choice of gratings and other requirements

