

# Argo-POWER<sup>SIM</sup> Slide V2

Fluorescent calibration slide with integrated power meter for high magnification systems, designed for the quality control of microscopes (QC).

**The SIM core of the POWER<sup>SIM</sup> slides are specifically designed designed for structured illumination microscopes.**

**Example of suited systems**

# 昊量光电



# 昊量光电

The Argo-POWER line of products integrates a compact yet sensitive power sensor. No need to modify your setup. Its shape and thinness make it fit adequately inside your microscope.

No more “Will 20% laser power kill my cells”. Using Daybook, our companion software, you can measure live power, record timelapses and trace input vs on-sample power graph. The software can generate ready-to-print reports based on your results.

Track the aging of your light sources with a reliable device. Argo-POWER sensors are calibrated and traceable to National French metrology institute (LNE) standards. Each comes with a certificate. Re-certification services are available.

#### Get quality and performance assessment solutions

Contact your local product specialist or someone from Argolight for more information.

[Set a meeting !](#)

# 昊量光电

## User stories

### Maria Smedh and Julia Fernandez-Rodriguez

*Respectively Researcher, and Head of the Centre for Cellular Imaging, Core Facility of the University of Gothenburg*

*“The intensity of the new slide is ~ 3-5 times brighter, depending on wavelength, compared to the 1st gen. This makes the patterns easier to find and it is easier to measure them in the orange-red and far-red wavelength regions. (...)The imaging direction is labeled on the new slide. This is nice, especially for unexperienced users, since it reduces mistakes.”*

Test carried out on Zeiss confocal (LSM 780) with Argo-SIM V2.

[Read the full testimony.](#)

## Dr. Patrick Then

*Microscopy specialist, Microverse Imaging Center, Friedrich Schiller University Jena*

*"Both on our wide-field and confocal microscopes, it was clear to see that the brightness of the slides has been nicely improved on the new generation of slides. (...) Being a new facility, we are currently setting up some testing protocols, which we intend to do regularly to make sure our microscopes have reliable performances (...) For that, the Argolight slides are a very nice tool, because the measures are reproducible and signal intensity doesn't change much over time. We're thus using the slide as a standard testing target at our facility."*

Test carried out on Zeiss Widefield (Elyra 7 Lattice SIM) and ZEISS Confocal (LSM 980) with Argo-SIM V2.

[Read the full testimony.](#)

## Jeroen Kole

*Product manager at Confocal.NL*

*"The RCM system is able to achieve 120 nm resolution which is hard to quantify using fluorescent beads. The Argo SIM V2 allowed us to visualize and quantify the resolution in a manner without introducing user-bias. We were clearly able to observe the 120 nm spacing between the lines in the image below. I noticed that the V2 version of the Argo SIM is much brighter than the first version. (...) for the first generation I used to push the laser to 50-60%, the increased brightness of the second-generation Argo SIM allowed me to obtain a similar SNR with only 5% laser power!"*

Test carried out on Re-scan Confocal Microscope 2 (RCM2) with Argo-SIM V2.


[Read the full testimony.](#)

昊量光电



## Argolight slides are designed to do routine quality assessments and reproducibility assessments of light microscopes.

They are made to improve reproducibility of light microscopy image data through quality control management of instruments (QC).





Argo-POWER<sup>SIM</sup> slides are re-usable long lasting fluorescent slides. While intensity may fluctuate with time, we warranty that they will be fluorescent for a lifetime




[How stable are Argolight Slides ?](#)

Argo-POWER<sup>SIM</sup> glass core fluorescence excitation spectrum is a continuum. All slides are compatible with any excitation from 250 nm to 650 nm. Emission is also a continuum starting about 10-20 nm after excitation and ending around 800 nm.



# Each Argo-POWER<sup>SIM</sup> glass core contains several fluorescent microscopic patterns

Each pattern is designed to be processed with the free companion software to detect aberrations and shift in microscopes performance. Users can catch issues before they can impact image quality and/or use the slides images to troubleshoot the image.





# Argo-POWER<sup>SIM</sup> specifications



昊量光电



Lifetime warrantied  
fluorescence presence

[more information  
on fluorescence stability](#)



**Excitation range:**  
continuum 250-650 nm



**Storage conditions:**  
room temperature (10-40 °C)  
and under normal relative humidity  
(20-70 % RH)



**Power measurement :**  
from 10 $\mu$ W to 100mW



**Timelapse compatible :**  
Measure several hours of data



**Dimensions:**

75x25x6 mm



**Emission range:**

continuum from the excitation  
wavelength plus 15 nm,  
to 800 nm



**Imaging technology**

**compatibility:**

any fluorescence-based imaging  
except depletion-based  
technology and  
multiphoton imaging

[more information on compatibility](#)



**Usable for wavelengths :**

from 350nm to 1100nm



**Materials:**

Anodized aluminum enclosure  
with an AG03 glass core



**Immersion medium compatibility:**  
dry, oil: no limitation,  
water objectives: less than 20 min at a time



**Light exposure damage threshold:**  
 $50 \text{ GW/cm}^2$  irradiance  
(peak or average)



Realtime power measurement


[More information in the user guide](#)

## Content of the slide

Each Argo-POWER<sup>SIM</sup> slide contains 27 fluorescent patterns.








### Target

PAT-AG03-EM2-A2

This pattern consists of concentric circles with increasing radii from 10  $\mu\text{m}$  to 120  $\mu\text{m}$  with a step of 10  $\mu\text{m}$ , featuring a target.



### Field of rings

PAT-AG03-EM2-B2

This pattern consists of a matrix of 21 $\times$ 21 rings, separated by 5  $\mu\text{m}$ , on a total field of 100  $\mu\text{m}$   $\times$  100  $\mu\text{m}$ . The field of rings is surrounded by eight landmarks and exhibits a 3  $\mu\text{m}$  long cross in its center.



### 4 $\times$ 4 Intensity gradation

PAT-AG03-EM2-C2

This pattern consists of sixteen 6  $\mu\text{m}$ -wide squares having different fluorescence intensity levels following a linear evolution, organized in a 4 $\times$ 4 matrix.



### 2 $\times$ 16 Intensity Gradation

PAT-AG03-EM2-D2

This pattern consists of twice sixteen 15.0  $\mu\text{m}$   $\times$  0.7  $\mu\text{m}$  rectangles having different fluorescence intensity levels following a linear evolution, organized in a 2 $\times$ 16 matrix.



### 3D Crossing stairs 1 $\mu\text{m}$ step

PAT-AG03-EM2-I5

This pattern consists of twice 11 empty cylinders embedded at different depths, like two crossing stairs, surrounded by four pillars. The step is: 1  $\mu\text{m}$ .



### 3D Crossing stairs 0.5 $\mu\text{m}$ step

PAT-AG03-EM2-I6

This pattern consists of twice 11 empty cylinders embedded at different depths, like two crossing stairs, surrounded by four pillars. The step is: 0.5  $\mu\text{m}$ .



### 3D Crossing stairs 0.25 $\mu\text{m}$ step

PAT-AG03-EM2-I7

This pattern consists of twice 11 empty cylinders embedded at different depths, like two crossing stairs, surrounded by four pillars. The step is: 0.25  $\mu\text{m}$ .



### 3D Crossing stairs 0.125 $\mu\text{m}$ step

PAT-AG03-EM2-I8


This pattern consists of twice 11 empty cylinders embedded at different depths, like two crossing stairs, surrounded by four pillars. The step is: 0.125  $\mu\text{m}$ .



### Word ARGOLIGHT

PAT-AG03-EM2-J2

This pattern consists of the letters forming the company name "Argolight", and surrounded by an  $80\ \mu\text{m} \times 18\ \mu\text{m}$  frame.



### 3D Matrix of rings

PAT-AG03-EM2-K1


This pattern consists in a 3D matrix of  $9 \times 9 \times 9$  rings, separated by  $5\ \mu\text{m}$ , on a total volume of  $40\ \mu\text{m} \times 40\ \mu\text{m} \times 40\ \mu\text{m}$ .



### Field of rings on a background

PAT-AG03-EM2-L1


This pattern consists in a matrix of  $9 \times 9$  rings, separated by  $5\ \mu\text{m}$ , on a total field of  $40\ \mu\text{m} \times 40\ \mu\text{m}$ , on a fluorescent background that is  $10\ \mu\text{m}$  below.



### Grid

PAT-AG03-EM2-N1


This pattern consists of a grid with a size of  $110\ \mu\text{m} \times 110\ \mu\text{m}$  and a step of  $10\ \mu\text{m}$ , containing crosses of  $5\ \mu\text{m}$  length in five squares around the center.



Geometrical figure:  
triangle

PAT-AG03-EM2-M2

One of the geometrical figures: a triangle.



Geometrical figure:  
circle

PAT-AG03-EM2-M1

One of the geometrical figures: a circle.



Geometrical figure:  
heptagon

PAT-AG03-EM2-M6

One of the geometrical figures: a heptagon.




Geometrical figure:  
hexagon

PAT-AG03-EM2-M5

One of the geometrical figures: a hexagon.






**Geometrical Figure**  
**octagon**

PAT-AG03-EM2-M6

One of the geometrical figures: an octagon.



**Geometrical figure**  
**pentagon**

PAT-AG03-EM2-M4

One of the geometrical figures: a pentagon.



**Geometrical figure:**  
**square**

PAT-AG03-EM2-M3


One of the geometrical figures: a square.



**Geometrical figure:**  
**star**

PAT-AG03-EM2-M8


One of the geometrical figures: a star with 16 arms.



### Gradually spaced lines

PAT-AG03-EM2-E8

This pattern consists of pairs of 36  $\mu\text{m}$ -long lines whose spacing gradually increases, from 0 nm to 390 nm, with a step of 30 nm. One set of lines is present: ascending (+45°).



### Gradually spaced lines

PAT-AG03-EM2-E5

This pattern consists of pairs of 36  $\mu\text{m}$ -long lines whose spacing gradually increases, from 0 nm to 390 nm, with a step of 30 nm. One set of lines is present: horizontal.



### Gradually spaced lines

PAT-AG03-EM2-E6


This pattern consists of pairs of 36  $\mu\text{m}$ -long lines whose spacing gradually increases, from 0 nm to 390 nm, with a step of 30 nm. One set of lines is present: vertical.



### Gradually spaced lines

PAT-AG03-EM2-E7

This pattern consists of pairs of 36  $\mu\text{m}$ -long lines whose spacing gradually increases, from 0 nm to 390 nm, with a step of 30 nm. One set of lines is present: descending (-45°).




### Matrix of crosses

PAT-AG03-EM2-F2

This pattern consists of a matrix of 4×4 crosses, having a length of 5  $\mu\text{m}$  and a step of 10  $\mu\text{m}$ , surrounded by a 40  $\mu\text{m}$ -wide frame.

The crosses are composed of vertical lines that are in the same plane, and by horizontal lines, going gradually deeper within the glass.


The spacing between the vertical and horizontal lines gradually increases, from 0.1  $\mu\text{m}$  to 1.6  $\mu\text{m}$ , with a step of 0.1  $\mu\text{m}$ .



### Sphere

PAT-AG03-EM2-G2

This pattern consists of three circles with a diameter of 25  $\mu\text{m}$  in different orthogonal planes, featuring the equator and two meridians of a sphere.



### Repositioning crosses

PAT-AG03-EM2-H2

The repositioning crosses are 20  $\mu\text{m}$  long.

[More information in the user guide](#)

Argo-POWER<sup>SIM</sup> Slide V2

## Find product documentation

Find the right documentation, user guide or information about your Argolight product.

[Knowledge Center](#)

昊量光电

[User guide](#)

[Software guide](#)